

# POWER UP

YOUR WHEELCHAIR

## Introducing the Alber **E-PILOT**

**The e-pilot turns your wheelchair into a sporty vehicle in a matter of seconds.**

In just a few simple steps, you can attach the e-pilot to your wheelchair to create a dynamic combination for everyday life. Because of its efficient Lithium-Ion battery and ergonomic design, this powerful add-on is an eco-friendly mobility alternative with maximum driving pleasure.

**Range and smoothness:**

Travel up to 31 miles on a single battery charge.

**Speed:**

Basic version 3.7 mph, optionally extendable to 12.4 mph through the Mobility Plus Package.

**More driving comfort:**

Stand with patented swivel function to navigate obstacles or uneven terrain.

**Mobility with style:**


Integrated battery pack and internal cable routing.

**Easy to handle:**

Easy operation using two control elements on the handles.

**Manufacturer:**

Alber





# POWER UP

YOUR WHEELCHAIR

## FRONTWHEEL

The RGK FrontWheel takes the sting out of pushing over rough terrain.



The RGK FrontWheel is a great add-on for the outdoorsy wheelchair user. It makes activities like walking the dog or camping a piece of cake. The way it attaches to an RGK wheelchair adds strength to the front end of the chair, whilst the front wheels are lifted off the ground. It has a much smaller turning circle than most of its competitors and is designed to minimise the impact on your chair's overall manoeuvrability.

**Weight:**  
2.6 kg

**Maximum user weight:**  
125 kg

**Solid and Stable:**  
200x50 pneumatic tyre with anti-flutter fork.

**Manufacturer:**  
RGK



## ATTITUDE hybrid

With 10 manual gears and an intuitive lightweight motor, the Hybrid is the preferred blend for many. It can be fitted to a number of different chairs and once it has been set up, you can effortlessly dock it onto your chair in about 6 seconds. The level of power assistance is adjustable, seamlessly kicking in at the resistance level you've set. This is a real game changer if you hit a steep hill or get caught out in the wrong gear. The built-in motor also helps to improve the front wheel's grip.

Experience freedom like never before.



**Battery:**  
14.5 Ah Lithium-Ion battery for maximum range performance.

**Motor:**  
The powerful 250 watt motor supports a top speed of up to 14.5 mph, ideal for hilly terrain or longer distances.

**Manufacturer:**  
Sunrise Medical



# LIPO LOMO micro

**A compact 12-inch  
cruising bike,  
perfect for travelling  
and cities.**



Stricker have built a handbike even more compact than the original Lipo Lomo. With the Lipo Lomo Micro you get a lightweight, but powerful add-on that will soon become an essential part of your everyday life. The Micro's whisper-quiet electric motor can take you up to 15 miles on a single charge, which makes it ideal for shopping trips and sightseeing tours.

**STRICKER**

**Range:**

The compact Lithium-Ion battery allows you to travel up to 15 miles when fully charged.

**Motor:**

The 250 watt 36 V Motor allows for a speed of nearly 7.5 mph.

**Weight:**

With battery 14.9 kg and without battery 13 kg.

**Manufacturer:**

Stricker

**STRICKER**



# E-FIX

**Give your existing lightweight  
wheelchair powerchair-like  
performance with the Alber e-fix.**



The Alber e-fix is a lightweight power add-on solution that can be fitted to most types of manual wheelchair. It allows the user to tackle hills and longer distances with ease. The built-in motors provide all the propulsion, so you can just sit back and steer using the ergonomically designed controller.



**e-fix®**



**Speed:**

0.37 to 3.7 mph

**Total weight:**

18.9 kg

**Max user weight:**

120 kg

**Manufacturer:**

Alber

**alber**



# POWER UP

YOUR WHEELCHAIR

## Introducing the *SmartDrive* **MX2+**



**The world of power assist products has been revolutionised over the past few years and Max Mobility is at the forefront of this revolution. Meet the SmartDrive MX2+: one of the most compact and portable power assist systems out there.**

The SmartDrive MX2+ can be fitted on both folding and rigid wheelchairs and enables you to cover longer distances and gradients with ease.

It also comes with a fashionable PushTracker, which is worn around the wrist and keeps track of such stats as push count and coast time. Using the integrated

PushTracker app, you can monitor your daily activities and also further customise this versatile power assist.

### **Range and smoothness:**

This small, yet substantial add-on can tackle 12 miles on a single charge.

### **Battery:**

A powerful Lithium-Ion (LiFePO4) battery supplies 36 V/3.4 Ah of power.

### **Weight:**

At a featherlight 5.7 kg, the SmartDrive MX2+ is one of the lightest power add-ons out there.

### **Effortless and versatile:**

The SmartDrive sits on the back of your chair and can be fit on nearly any make or type of wheelchair.

### **Manufacturer:**

Max Mobility

*SmartDrive*  
now part of permobil

