

Avantgarde XXL 2



Sporting challenges or active everyday life – with the Avantgarde XXL 2, you can cross all boundaries. The active wheelchair offers the ideal combination of mobility and stability.

- For bariatric patients with various indications; many options available
- Particularly stable
- Max. load capacity of 180 kg
- Seat upholstery with reinforced insert
- Side panels with neoprene insert
- Single-panel footplate, angle-adjustable, swing-away and removable
- 24" or 26" rear wheels with quick-release axle, wider hub, reinforced hollow rims and aluminium push rings
- Rear wheel adapter 0°
- Low net weight
- Aluminium tube frame

Seat width	46 – 62 cm
Seat depth	40 – 56 cm
Seat height front/rear	46 – 53/41 – 53 cm
Width	66.5 – 82 cm
Length	79 – 114 cm
Overall height	70 – 103.5 cm
Folding size from	32 cm
Lower leg length	32 – 49 cm
Variable armrest height	20 – 30 cm
Max. user weight	180 kg
Weight approx.	17 kg